



▶ Evaporating dishes

1 | Nahita evaporating dishes are ideal for drying solids in stoves or directly to the flame of a Bunsen burner, since they are made of porcelain resistant to extreme temperatures up to 1050° C.

2 | They are available in two different models and with several capacities and are glazed inside

▶ Round bottom evaporating dishes




1 | They present spout and are glazed inside and partially outside.

Code	Capacity	Diameter	Height
CNG001	15 mL	47 mm	20 mm
CNG002	35 mL	60 mm	25 mm
CNG003	50 mL	70 mm	28 mm
CNG004	60 mL	74 mm	30 mm
CNG005	75 mL	82 mm	32 mm
CNG006	100 mL	90 mm	35 mm
CNG007	125 mL	98 mm	38 mm
CNG008	140 mL	100 mm	39 mm
CNG009	150 mL	104 mm	40 mm
CNG010	200 mL	110 mm	44 mm
CNG011	270 mL	120 mm	49 mm
CNG012	300 mL	126 mm	52 mm
CNG013	400 mL	150 mm	58 mm
CNG014	1000 mL	206 mm	64 mm

▶ Flat bottom evaporating dishes (casseroles)



DIMENSIONS - Flat bottom



1 | They present spout and handle, and are glazed inside and outside except for rim and outside bottom.

Code	Capacity	Diameter	Height
CNK001	30 mL	50 mm	30 mm
CNK002	70 mL	66 mm	39 mm
CNK003	150 mL	90 mm	48 mm
CNK004	200 mL	96 mm	52 mm
CNK005	250 mL	108 mm	60 mm
CNK006	500 mL	132 mm	70 mm