



► Porcelain crucibles

1 | Nahita offers this wide range of porcelain crucibles with high, medium or tall form that are ideal for incineration or calcination of substances in analytical chemistry laboratories


2 | Crucibles are totally glazed except for rim and outside bottom and can withstand extreme temperatures up to 1050° C without cracking or deforming

3 | Each crucible is supplied complete with a glazed lid that presents a handle for an easy manipulation and that allows the exit of vapours or gases produced during the calcination reaction.

► High form crucible



DIMENSIONS - High form crucible




1 | With lid. Glazed interior and exterior except rim and base.

Code	Capacity	Diameter	Height
CPB001	15 mL	32 mm	34 mm
CPB002	18 mL	34 mm	36 mm
CPB004	20 mL	36 mm	38 mm
CPB005	25 mL	38 mm	40 mm
CPB006	30 mL	40 mm	46 mm
CPB007	40 mL	42 mm	54 mm
CPB008	50 mL	43 mm	58 mm

► Medium form crucible



DIMENSIONS - Medium form crucible




1 | With lid. Glazed interior and exterior except rim and base.

Code	Capacity	Diameter	Height
CPC001	5 mL	25 mm	21 mm
CPC002	10 mL	30 mm	28 mm
CPC003	15 mL	35 mm	30 mm
CPC004	25 mL	40 mm	36 mm
CPC005	30 mL	42 mm	41 mm
CPC006	40 mL	48 mm	42 mm
CPC007	50 mL	53 mm	46 mm
CPC008	70 mL	56 mm	56 mm
CPC009	100 mL	63 mm	59 mm
CPC010	300 mL	90 mm	78 mm

► Low form crucible



DIMENSIONS - Low form crucible



1 | With lid. Glazed interior and exterior except rim and base.

Code	Capacity	Diameter	Height
CPD001	15 mL	43 mm	23 mm
CPD002	25 mL	47 mm	27 mm
CPD003	30 mL	51 mm	30 mm
CPD004	50 mL	59 mm	38 mm