



Desiccators are used for the desiccation of moist substances or for the storage of materials sensitive to air humidity. Their construction enables usage of vacuum processes. Both, the contact surface of the lid and the bottom of desiccator are precision ground. All parts of the desiccator are exchangeable and therefore can be provided as

replacement parts. The different design of particular parts enables the user to compile desiccator as necessary. Please pay close attention to specified diameter of combined parts.


► Desiccators with plastic knob



The desiccators are made in borosilicate glass perfectly transparent. Lid with plastic knob included. Porcelain disk not included.

Code	Inner diameter	External diameter
BDD006	100 mm	151 mm
BDD007	150 mm	210 mm
BDD008	200 mm	269 mm
BDD009	250 mm	329 mm
BDD010	300 mm	392 mm

dimensions



► Desiccators with nz 24/29 and valve



Desiccator is made from Simax 3.3 borosilicate Glass. It is designed for a variety of vacuum work. All parts are exchangeable and can be provided as replacement parts.

Typical use: desiccation of samples, storage of materials sensitive to humidity

Code	d <sub>1</sub>	d <sub>2</sub>	h
BDB020	153 mm	94 mm	177 mm
BDB021	215 mm	143 mm	247 mm
BDB022	270 mm	192 mm	300 mm
BDB023	320 mm	239 mm	349 mm
BDB024	380 mm	285 mm	424 mm



Code	d <sub>2</sub>	d <sub>1</sub>
BDB011	100 mm	151 mm
BDB012	150 mm	210 mm
BDB013	200 mm	269 mm
BDB014*	250 mm	329 mm
BDB015	300 mm	392 mm

\* este artículo se suministra bajo pedido




▶ Desiccators with glass knob



ISO  
13130

dimensions



Desiccator is made from Simax 3.3 borosilicate Glass. It is designed for a variety of vacuum work. All parts are exchangeable and can be provided as replacement parts.

Typical use: desiccation of samples, storage of materials sensitive to humidity.

Code	$d_1$	$d_2$	h
BDD011	153 mm	94 mm	191 mm
BDD012	215 mm	143mm	261 mm
BDD013*	270 mm	192 mm	311 mm
BDD014	320 mm	239 mm	363 mm
BDD015	380 mm	285 mm	438 mm

\* these references are available under request


▶ Desiccator with straight stopcock



ISO  
13130

BORO  
3.3

dimensions



The desiccators are made in borosilicate glass perfectly transparent. Lid with glass stopcock included. Porcelain disk not included.

Code	$d_2$	$d_1$
BDB006	100 mm	151 mm
BDB007	150 mm	210 mm
BDB008*	200 mm	269 mm
BDB009	250 mm	329 mm
BDB010	300 mm	392 mm

\* these references are available under request

▶ Bell jar with knob



BORO  
3.3

Code	Ø base	h
BDJ001	150 mm	150 mm
BDJ002*	150 mm	250 mm
BDJ003	150 mm	200 mm
BDJ004	200 mm	200 mm
BDJ005	200 mm	300 mm
BDJ006	300 mm	400 mm

\* este artículo se suministra bajo pedido



▶ ACCESSORIES FOR DESSICATORS

▶ PP knob with hole



Code	Knob diameter	Hole diameter
BDL028	52 mm	9 mm

▶ PP knob for knob lid



Code	Base diameter
BDL046	52 mm

▶ Lids for desiccators with nz 24/29 for knob lids



BORO  
3.3

Code	d	h	Ground joint
BDL018	153 mm	64 mm	24/29
BDL019	215 mm	89 mm	24/29
BDL020	270 mm	97 mm	24/29
BDL021	320 mm	111 mm	24/29
BDL022	380 mm	136 mm	24/29

▶ Lids for desiccators with glass knob



BORO  
3.3

Code	d	h
BDL006	153 mm	78 mm
BDL007	215 mm	103 mm
BDL008	270 mm	108 mm
BDL009	320 mm	125 mm
BDL010	380 mm	150 mm



BORO  
3.3

Code	For desiccator	Base diameter	h
BDL036	100 mm	153 mm	83 mm
BDL037	150 mm	213 mm	95 mm
BDL038	200 mm	272 mm	123 mm
BDL039	250 mm	333 mm	139 mm
BDL040	300 mm	395 mm	140 mm



BORO  
3.3

Code	For desiccator	Base diameter	h
BDL041	100 mm	153 mm	84 mm
BDL042	150 mm	213 mm	96 mm
BDL043	200 mm	272 mm	124 mm
BDL044	250 mm	333 mm	140 mm
BDL045	300 mm	395 mm	141 mm



▶ ACCESORIES FOR DESSICATORS

▶ Lids for desiccators with straight stopcock



dimensions

$d_1$ : base diameter  
 $d_2$ : hole inner diameter  
 $h$ : height



Code	For desiccator	$d_1$	$d_2$	$h$
BDL029*	100 mm	153 mm	9 mm	84 mm
BDL030*	150 mm	213 mm	9 mm	96 mm
BDL031*	200 mm	272 mm	9 mm	124 mm
BDL032	250 mm	333 mm	9 mm	140 mm
BDL033	300 mm	395 mm	9 mm	141 mm

\* este artículo se suministra bajo pedido

▶ Stopcock 24/29 for desiccator



Code	Ground joint
BDL023	24/29

▶ Straight desiccator stopcock with joint



Code	Base diameter
BDL034	8 mm

▶ Stopcock 24/29 for dessicator



Code	Ground joint
BDL035	24/29