



BJJ - For tubes, PVC

1 | Nahita offers a wide variety of different size and capacity test tube PVC racks.


2 | Racks are an essential element in any laboratory to keep samples perfectly tidy and hold them in fridges, freezers or on the bench top.

▶ For tubes of 13 mm Ø



Code	Alignment	N° tubes	Dimensions (LxHxW)
BJJ001	2x3	6 tubes	55x45x40 mm
BJJ002	1x6	6 tubes	112x38x40 mm
BJJ003	2x6	12 tubes	118x48x40 mm
BJJ004	1x12	12 tubes	250x44x40 mm
BJJ005	2x12	24 tubes	255x48x40 mm
BJJ006	1x24	24 tubes	490x50x40 mm
BJJ007	3x16	48 tubes	310x75x40 mm
BJJ009	5x12	60 tubes	200x131x40 mm
BJJ010	10x10	100 tubes	195x195x40 mm
BJJ011	7x20	140 tubes	410x150x40 mm
BJJ012	10x15	150 tubes	291x195x40 mm

▶ For tubes of 16 mm Ø



Code	Alignment	N° tubes	Dimensions (LxHxW)
BJJ013	2x6	12 tubes	160x60x40 mm
BJJ014	2x12	24 tubes	268x65x40 mm
BJJ015	2x6	12 tubes	160x60x65 mm
BJJ016	2x12	24 tubes	268x65x65 mm

▶ For tubes of 22 mm Ø



Code	Alignment	N° tubes	Dimensions (LxHxW)
BJJ017	2x6	12 tubes	182x68x65 mm
BJJ018	2x12	24 tubes	350x82x658 mm
BJJ019	2x3	6 tubes	185x65x90 mm
BJJ020	2x6	12 tubes	182x68x90 mm
BJJ022	2x12	24 tubes	350x83x90 mm



BJJ - For tubes, PVC

▶ For tubes of 32 mm Ø



Code	Alignment	N° tubes	Dimensions (LxHxW)
BJJ023	2x3	6 tubes	130x90x65 mm
BJJ025	2x4	8 tubes	185x95x65 mm
BJJ026	2x6	12 tubes	250x85x65 mm

▶ For tubes of 40, 50 and 60 mm Ø



Code	Alignment	N° tubes	Ø tubes	Dimensions (LxHxW)
BJJ027	2x3	6 tubes	40 mm	177x120x60 mm
BJJ028	2x4	8 tubes	40 mm	250x130x60 mm
BJJ029	2x2	4 tubes	50 mm	150x150x60 mm
BJJ030	2x3	6 tubes	50 mm	220x150x60 mm