



Porcelain spot plates - Nahita spot plates are made of porcelain and are highly resistant to acids and bases and thermal shock.



1| Thanks to their great chemical resistance are ideal to carry out small chemical reactions while due to their white colour are also perfect for colour reactions or even for staining small specimens.

2| Plate surface and sides are glazed to ease cleaning and each well is approximately 20 mm diameter and 5 mm depth.

Code	N° wells	Lenght	Widht
BPL014	6	85 mm	54 mm
BPL015	12	115 mm	90 mm