



► Irrigator barrel shaped



- 1 | Tubulated near the bottom with loose drop on cover with knob.
- 2 | Provided with standardized ground neck near bottom and bleeding cock.

Code	Ground	Capacity	ø lid
CLB001	29/32	3L	100 mm
CLB002	29/32	5 L	100 mm
CLB003	29/32	10 L	127 mm
CLB004	29/32	20 L	160 mm

Accessories



CODE: CLB005



CODE: CLB006

Code	Description
CLB005	RIN 29/32 Stopcock for barrel
CLB006	Lid for 5L barrels