




▶ Sponge scrubber Mousse range N°1



- 1 | Autoclavable
- 2 | Made from fibres
- 3 | For all kind of test tubes with diameter < 25 mm

Code	Sponge width	Sponge length	Total length	Application
ZCF003	22 mm	170 mm	220 mm	Test tubes Ø < 25 mm

▶ Sponge scrubber Mousse range N°2, 3 and 4



- 1 | Made from polyethylene foam that is durable and long lasting
- 2 | 20 times more resistant than a standard sponge or brush
- 3 | Do not absorb water
- 4 | Autoclavable. Easily cleaned with hot or cold water

Code	N°	Sponge width	Sponge length	Total length	Application
ZCF002	2	70 mm	100 mm	250 mm	Beakers, flasks, jars, etc.
ZCF001	3	70 mm	100 mm	480 mm	Measuring cylinders, etc.
ZCF004	4	20 mm	380 mm	440 mm	Test tubes, Erlenmeyer flasks, separating funnels, etc.,

▶ Flexible brushes



- 1 | Flexible rubber hose covered with a sheath of polyvinyl chloride

Code	Ø bristle	Bristle length	Total length	Material
ZCF005	30 mm	160 mm	450 mm	Nylon
ZCF006	45 mm	200 mm	200 mm	Nylon