




► Shoe covers

- 1 | PE plastic shoe covers. Ambidextrous.
- 2 | With elastic at ankle height. White colour.
- 3 | For general use against minor risks. Single use.
- 4 | 100% Latex free product.

Protective gloves for general use, especially suitable for use in the Healthcare, Medical, Hospital, Laboratory Equipment, Chemical, Food, Painting, Industrial, Automotive, Construction and Agriculture sectors.

CE-certified:

- EU 2016/425: Relating to personal protective equipment. PPE Category I.
- EU 2017/745: For Medical Devices, Class I.
- EU 10/2011: On plastic materials intended to come into contact with foodstuffs, considering their migration limits.
- EC 1935/2004: Art. 5, paragraph I.



Code	Size	Dimensions	Thickness	Colour
RBN001	U	410 x 155 mm	0,018 mm	White
RBN002	U	410 x 155 mm	0,018 mm	Green
RBN003	U	410 x 155 mm	0,018 mm	Blue