



► Nitrile gloves

- 1 | 100% unpowdered Nitrile glove. Highly elastic. Ambidextrous.
- 2 | For general use against minor hazards. Single use.
- 3 | Micro-textured outer surface on fingertips for better grip. -Smooth inner surface to facilitate grip and reduce dermatitis.
- 4 | Reinforced cuff.
- 5 | Phthalate-free product.

Protective gloves for general use, especially suitable for use in the Healthcare, Medical, Hospital, Laboratory Equipment, Chemical, Food, Painting, Industrial, Automotive, Construction and Agriculture sectors.

CE certified:

- MDR 2017/745: Class I. Medical Devices
- EU 2016/425: Category III. Related to Personal Protective Equipment.

Regulations concerning glove protection:

- EN-420:2003: General protection
- EN-374 (374.1, 374.2, ISO 16523-1:2015, 374.4, 374.5): Against chemicals and micro-organisms. Permeability.
- EN-455 (455.1, 455.2, 455.3, 455.4): Protective gloves for exploration.



Code	Size	Weight	Colour
RBL005	S 6/7	5	Cobalt Blue
RBL006	M 7/8	5	Cobalt Blue
RBL007	L 8/9	5	Cobalt Blue
RBL008	XL 9/10	5	Cobalt Blue