




▶ Hooded coverall

Protective cover for general use, especially suitable for use in the Healthcare, Hospital, Laboratory Equipment, Chemical, Aesthetics, Food, Cleaning, Painting, Industry and Agriculture sectors.

- Single use long sleeve coverall made of Polypropylene,
- Zip fastening.
- Elasticated hood and cuffs.
- Hood for head protection.
- Individually bagged.
- 100% Latex free product.

CE certified

1. EU 2016/425: Relating to personal protective equipment. PPE Category I.
2. EU 2017/745: For Medical Devices, Class I.
3. EU 10/2011: On plastic materials intended to come into contact with foodstuffs, considering their migration limits.
4. EC 1935/2004: Art. 5, paragraph I.



Code	Size	Colour
RBC001	M	White
RBC002	L	White
RBC003	XL	White
RBC004	XXL	White
RBC005	XXXL	White