




► Polyethylene Apron

- Polyethylene (LDPE) apron, with head opening and fastening straps on the back.
- For general purpose protection against minor hazards. Single use.
- 100% Latex free product.

CE-certified:

- EU 2016/425: Relating to personal protective equipment. PPE Category Category I.
- EU 10/2011: On plastic materials intended to come into contact with foodstuffs, studying their migration limits.
- EC 1935/2004: Art. 5, section I.



| | |
|-------------------|-------------|
| Code | RBK001 |
| Size | U |
| Dimensions | 120 X 78 cm |
| Thickness | 0,02 mm |
| Colour | White |