



▶ Graduated cone shape thick glass



- 1 | Soda glass
- 2 | Graduated cones are conical, present graduation.
- 2 | They are used for measuring and transferring approximate volumes of liquids.
- 3 | The base is wide and stable and they are also provided with a spout for an easy pouring.
- 4 | Graduated cones are available in different capacities to adapt to users' necessities.

Code	Height	Lower diameter	Upper diameter	Capacity
CBD004	103 mm	22 mm	50 mm	60 mL
CBD005	119 mm	34 mm	80 mm	125 mL
CBD006	144 mm	40 mm	93 mm	250 mL
CBD007	189 mm	46 mm	120 mm	500 mL

GRADUATED CONE SHAPE THICK GLASS  
[C] BEAKERS, FLASKS AND OTHER CONTAINERS

Note: the manufacturer reserves the right to modify the characteristics and photographs of the products