



## Macropipette tips

Code MJQ021 - MJQ022

- Autoclavable, made of Polypropylene.
- Productions in controlled atmosphere controlada, with personnel protected with epis, gown, gloves, mask. The measurements may differ somewhat due to the continuous molding process.
- Hazards data:
  - Health: 0 Not dangerous
  - Flammable: 1 Flammable (Class A)
  - Reactivity: 0 Non-reactive



## Technical specifications

Reference	MJQ021	MJQ022
Type	A	B
Compatible with pipette	MGB008 (Digipette Standard)	MGB095 (Digipette Premium)
Maximum volumen		5000 $\mu$ L
Graduation	1000 $\mu$ L, 2000 $\mu$ L, 3000 $\mu$ L, 4000 $\mu$ L y 5000 $\mu$ L	

## Technical drawing

