

Issuer's Name / Producer : GLASSCO LABORATORY EQUIPMENTS PVT. LTD.

Issuer's Address / Producer : Manglai, P.O. Khudda Kalan, Ambala Cantt - 133001
Haryana, India


Object of Declaration : Pear shape flask

Material : Borosilicate Glass

Purpose of Use : Pear-shaped flasks are used for evaporating solutions to dryness post synthesis using a rotary evaporator, the 'rounded V' shape of the flasks enables solid materials to be scraped out more Efficiently than from a round bottom flask.

Specification of Flask specified in Table-1, DIN 12383

Catalogue No.	Product Description	Socket Size	Capacity ml	Pack Qty
PBM001	Pear shape flask	14/23	25	10
PBM005	Pear shape flask	14/23	50	10
PBM006	Pear shape flask	14/23	100	10



*non DIN

The object of declaration described above is in conformity with the requirements of the following documents:

- **DIN ISO 3585** Borosilicate Glass 3.3
- **DIN 12383** Laboratory Glassware - Pear shape Flask

Material: The above are made from DIN ISO 3585:1998, heat resistant Boro 3.3 Glass that meet ASTM E-438, TYPE-1, CLASS- A Standard.

Compliance: These flask comply in every respect characteristics, measures and accuracy are in full Conformity With the provisions of the standard DIN 12383.

Additional Information:

Glassco has installed the quality control system as per ISO 9001 and products delivered to the market are in full conformity with the technical Specification as per standard and with all fundamental requirements to such products.

Certificate No. IPK 0039.69-01 issued by intertek certification limited, A UKAS Management System Accredited Body

