



▶ Pocket balance, 5043 series

Nahita pocket scales have a modern and very compact design, making them easy to carry in a handbag or even in your pocket. They allow precise weighing of for example precious stones and metals, medicinal herbs or ingredients for a diet:

1 | Operation via backlit touch screen for easy reading.

- 2 | Compact and rugged design with plastic cover to protect the sensor from shocks and overloads, and stainless steel weighing pan.
- 3 | Weighing, piece counting and tare function via hotkey.
- 4 | One-point calibration (calibration weight not included).
- 5 | Automatic shut-off after 1 minute of non-use to extend battery life.
- 6 | Low battery (LO) and overload (O_Ld) indicators on display.



Code	KBH001	KBH002	KBH003
Capacity	100 g	200 g	500 g
Sensibility	0.01 g	0.01 g	0.01 g
Measuring units	g, ct, dwt, oz, ozt		
Functions	Weighing, piece counting and tare functions.		
Plate dimensions	60x60 mm		
Dimensions	70x11 x 0 mm		
Power	2 batteries 1,5V Type: AAA		