



Pipeteur électronique

RÉF. - MBD002

- Intuitive and convenient speed with 6 levels
- Lightweight pipette filler compatible with glass or plastic serological pipettes (0.1-100mL)
- Ergonomic design for fatigue-free and effective streamlined pipetting
- LCD screen displays battery status and pipetting speed
- Quiet and powerful motor and pup fills 25mL pipette in less than 5 seconds
- Long-lasting lithium-ion battery



SPÉCIFICATIONS TECHNIQUES

Référence	MBD002
Vitesse d'aspiration	25mL < 5s (6 niveau)
Vitesse de distribution	Moteur (6 changement) / La gravité
Batterie	Lithium-ion
Durée de vie de la batterie	Plus de 8 heures d'utilisation intermittente
Temps de charge	2-3 heures
Types de pipettes	Pipette en verre ou en plastique (0,1 à 100 ml), pipettes Pasteur
Filtre	0.45 µm Hydrophobe
Poids	200 g