





Combustion crucible - low form  
Creuset combustion, forme basse  
Crisol de combustión, forma baja



Code	D mm	H mm	V/ml approx.	Code JIPO	DIN
CPD005	30	19	5	641 321 100 301	DIN 12 904
CPD006	35	22	10	641 321 100 351	DIN 12 904
CPD007	40	25	17	641 321 100 401	DIN 12 904
CPD008	45	28	21	641 321 100 451	DIN 12 904
CPD009	50	32	34	641 321 100 501	DIN 12 904
CPD010	60	38	62	641 321 100 601	DIN 12 904
CPD011	70	44	91	641 321 100 701	DIN 12 904




Combustion crucible - middle  
Creuset combustion, forme moyenne  
Crisol de combustión, forma media



Code	D mm	H mm	V/ml approx.	Code JIPO	DIN
CPC011	30	25	10	641 321 200 301	
CPC012	35	28	12	641 321 200 351	DIN 12 904
CPC013	40	32	20	641 321 200 401	DIN 12 904
CPC014	45	36	30	641 321 200 451	DIN 12 904
CPC015	50	40	45	641 321 200 501	DIN 12 904
CPC016	60	48	80	641 321 200 601	DIN 12 904
CPC017	70	56	120	641 321 200 701	DIN 12 904
CPC018	80	62	200	641 321 200 801	




Combustion crucible - high form  
Creuset combustion, forme haute  
Crisol de combustión, forma alta



Sign	D mm	H mm	V/ml approx.	Code JIPO	DIN
CPB010	30	38	15	641 321 300 301	DIN 12 904
CPB011	35	44	26	641 321 300 351	DIN 12 904
CPB012	40	50	35	641 321 300 401	DIN 12 904
CPB013	45	56	50	641 321 300 451	DIN 12 904
CPB014	50	62	72	641 321 300 501	DIN 12 904
CPB015	60	75	120	641 321 300 601	DIN 12 904



Lid for combustion crucible  
Couvercle p/creuset  
Tapa para crisoles



Označení Sign	D 1 mm	D 2 mm	V/ml approx.	Číselník Code JIPO	DIN
CPF003	30	34		641 321 000 301	DIN 12 904
CPF004	35	39		641 321 000 351	DIN 12 904
CPF005	40	44		641 321 000 401	DIN 12 904
CPF006	45	49		641 321 000 451	DIN 12 904
CPF007	50	54		641 321 000 501	DIN 12 904
CPF008	60	64		641 321 000 601	DIN 12 904
CPF009	70	74		641 321 000 701	DIN 12 904
CPF010	80	84		641 321 000 801	