



ROGO-SAMPAIC
matériel de laboratoire



Water Distillation Unit BASIC PH4 4 lit / h

Code : 57 000 003

Technical features :

BASIC PH4 model is designed for reliable continuous operation incorporating a host of features unmatched by comparable stills. It is easy to use, safe to operate, thoroughly reliable and above all Low Cost. Heater is a 3KW chromium plated heater fitted with a safety thermostat. The boiler is made of high quality borosilicate glass.

All units carry CE marking.

Conductivity (using raw water) : 3.0 - 4.0 $\mu\text{s}/\text{cm}$

Conductivité (using treated water) : 1.5 - 2.0 $\mu\text{s}/\text{cm}$

pH : 5.5 - 6.0

Dimensions : 50 x 15 x 45 cm

Accessoires de rechange :

1	Heater	ref. : 57 000 100
2	Boiler	ref. : 57 000 104
3	Condenser	ref. : 57 000 105
4	Tubing set	ref. : 57 000 103
5	Joint	ref. : 57 000 101
6	Metal stand	ref. : 57 000 107

Distributed by :