

## COMECTA Chairs for laboratory

NEW

### FEATURES

Chairs for laboratory with adjustable height of 60 to 80cm. Seat made of polyurethane foam. Available with or without wheels and in aluminium or chromed.

### MODELS

Part No.	Ø base cm	Made in	Wheels	Adaptable weight cm
5609601	40	Aluminium	Yes	60 to 80
5609602	40	Chromed	Yes	60 to 80
5609603	40	Aluminium	No	60 to 80
5609604	40	Chromed	No	60 to 80



## COMECTA Transport Refrigerator HJ-1

DESIGNED TO TRANSPORT AND PRESERVE LABORATORY BIOLOGICAL SAMPLES.

MEDICAL DIRECTIVE CERTIFICATE: CLASS 1-93/43 EEC

### FEATURES

Strap and handles for transportation.

Opening on both sides.

Supplied with 10 packs of cooling gel. Easy access by opening at the top. Extra aluminium insulation.

Drainage.

Integrated thermometer with LCD screen for reading the temperature.



NEW

### MODEL

Part No.	Height / Width / Depth Interior cm			Height / Width / Depth Exterior cm			Capacity litres without gel	Capacidad litres with gel	Weight without gel Kg	Weight with gel Kg
5601057	25	33	17	32	39	27	14	4,2	4,5	9,8