

Stomacher® 400 Circulator Bags

The **Stomacher® 400 Bag** Range caters for both the Circulator bag (round-bottom) designed blender bags, and the Classic (square-bottomed) for the older square paddle designed Stomachers®.

Our blender bags are made from high quality food grade virgin LDPE granule, extruded to 65µm film thickness. All Stomacher® Bags are irradiated sterile. A certificate of irradiation is available upon request.



Standard bags
177 x 305mm
80-400ml
500 (sachets of 50)

ZHD508



Closure Bag
177mm x 305mm
80-400ml
250 (sachets of 10)

ZHD510



Strainer Bag
177mm x 305mm
80-400ml
200 (sachets of 50)

ZHD5009