



U-tubes

Code **BGJ001, BGJ008, BGJ009, BGJ010, BGJ011, BGJ012, BGJ006, BGJ007**

U-shaped tubes are used to determine the pressure at which the gas contained in a vessel is at. For this purpose, the tube is filled with a high-density liquid, such as Hg, while one end of the tube is connected to the vessel containing the gas. The higher the gas pressure, the greater the difference in height between the two Hg branches. This height difference will determine the pressure at which the gas is found (in mm Hg).

Endo Glassware U-tubes are made of borosilicate glass.



Code	Length (mm)	Ø (mm)	Tube spacing (mm)	Rim
BGJ001	150	16	50	YES
BGJ008	140	14	60	YES
BGJ009	160	16	66	YES
BGJ010	180	18	66	YES
BGJ011	200	20	66 / 72	YES
BGJ012	200	20	75	YES
BGJ006	200	20	54	NO
BGJ007	200	24	80	NO