




Mortar unglazed
Mortero sin vitrificar
Mortier biscuit,



Code	D mm	H mm	V/ml approx.	Code JIPO	DIN
BBJ001	51	32	25	641 331 211 000	
BBJ002	56	36	30	641 331 211 100	
BBJ018	63	41	70	641 331 211 110	DIN 12 906
BBJ003	76	45	75	641 331 211 120	
BBJ004	83	50	110	641 331 211 200	
BBJ005	90	56	160	641 331 211 220	DIN 12 906
BBJ006	105	64	220	641 331 211 300	
BBJ007	125	71	400	641 331 211 320	DIN 12 906
BBJ019	139	82	500	641 331 211 400	
BBJ009	150	90	650	641 331 211 500	DIN 12 906
BBJ010	180	104	1000	641 331 211 600	DIN 12 906
BBJ011	224	127	2250	641 331 211 700	
BBJ012	268	160	4300	641 331 211 800	



Pestle unglazed
Pilon pour mortier standard biscuit,
Mano para mortero, cabeza sin vitrificar



Code	D mm	L mm	Code JIPO	DIN
BBN007	14	54	641 331 213 000	
BBN008	24	115	641 331 213 100	DIN 12 906
BBN009	30	135	641 331 213 200	DIN 12 906
BBN010	36	150	641 331 213 300	DIN 12 906
BBN011	42	175	641 331 213 400	DIN 12 906
BBN012	55	210	641 331 213 500	DIN 12 906
BBN013	74	250	641 331 213 600	DIN 12 906
BBN014	85	300	641 331 213 700	DIN 12 906
BBN015	101	370	641 331 213 800	




Pestle unglazed
Pistill unglasiert
Mano para mortero, cabeza sin vitrificar



Sign	D mm	L mm	Code JIPO
BBM001	18	68	641 331 212 100
BBM002	20	88	641 331 212 200
BBM003	25	96	641 331 212 300
BBM004	34	140	641 331 212 400
BBM005	37	158	641 331 212 500



Mortar glazed
Mortero vitrificado
K mprȳcp ĉ k _gjjĉ*



Code	D mm	H mm	V/ml approx.	Code JIPO	DIN
BBH001	51	32	25	641 331 211 001	
BBH002	56	36	30	641 331 211 101	
BBH003	63	41	70	641 331 211 111	DIN 12 906
BBH004	76	45	75	641 331 211 121	
BBH005	83	50	110	641 331 211 201	
BBH006	90	56	160	641 331 211 221	DIN 12 906
BBH007	105	64	220	641 331 211 301	
BBH008	125	71	400	641 331 211 321	DIN 12 906
BBH009	139	82	500	641 331 211 401	
BBH010	150	90	650	641 331 211 501	DIN 12 906
BBH011	180	104	1000	641 331 211 601	DIN 12 906
BBH012	224	127	2250	641 331 211 701	
BBH013	268	160	4300	641 331 211 801	




Pestle glazed
Pilon, émaillé
Mano para mortero, cabeza vidriada



Code	D mm	L mm	Code JIPO	DIN
BBN024	14	54	641 331 213 001	
BBN016	24	115	641 331 213 101	DIN 12 906
BBN017	30	135	641 331 213 201	DIN 12 906
BBN018	36	150	641 331 213 301	DIN 12 906
BBN019	42	175	641 331 213 401	DIN 12 906
BBN020	55	210	641 331 213 501	DIN 12 906
BBN021	74	250	641 331 213 601	DIN 12 906
BBN022	85	300	641 331 213 701	DIN 12 906
BBN023	101	370	641 331 213 801	



Mortar / grinding surface unglazed (BBJ) ; glazed (BBH)
Mortier / Biscuit (BBJ) ; émaillé (BBH)
Mortero/ sin vitrificar (BBJ); vidriado (BBH)



Code	D mm	H mm	V/ml approx.	Code JIPO
BBJ013	85	65	150	641 331 829 100
BBH014	85	65	150	641 331 829 101
BBJ014	103	75	300	641 331 829 200
BBH015	103	75	300	641 331 829 201
BBJ015	121	88	420	641 331 829 300
BBH016	121	88	420	641 331 829 301
BBJ015	152	95	730	641 331 829 400
BBH017	152	95	730	641 331 829 401
BBJ017	178	110	1200	641 331 829 500
BBH018	178	110	1200	641 331 829 501