



Density petroleum hydrometer, L-50 Series

Code **KLG067, KLG068, KLG069, KLG070, KLG071, KLG072, KLG073**

- Manufactured according to BS 718
- Reference temperature: 15°C
- Amplitude of the measuring range: 0.050 g/mL
- Scale division: 0.0005 g/mL
- Total length: 335 mm
- Factory calibrated
- The instrument is supplied without certificate



Technical specifications

Code	KLG067	KLG068	KLG069	KLG070	KLG071	KLG072	KLG073
Range (g/mL)	0.650-0.700	0.700-0.750	0.750-0.800	0.800-0.850	0.850-0.900	0.900-0.950	0.950-1.000