



Culture tubes

Code **BGG013**, **BGG016**, **BGG006**, **BGG015**

Round bottom, without stopper

Code **BGG013**, **BGG016**



- Made of soda glass.
- With screw neck and round bottom.

Code	Ø (mm)	Length (mm)	Thickness (mm)	Capacity (mL)
BGG013	16	160	1.00	22
BGG016	19.25	200	1.00	42

Round bottom, with aluminium cap

Code **BGG006**, **BGG015**



- Made of soda glass.
- Tube with screw neck and round bottom.
- Aluminium cap with rubber gasket.

Code	Ø tube (mm)	Tube length (mm)	Thickness (mm)	Capacity (mL)
BGG006	16	160	1.00	22
BGG015	19.25	200	1.00	42