




Staining dishes and baskets

For Coplin staining, Hellendahl staining, Schieferdecker staining and Borrel staining

Endo glassware is provided with a wide range of staining jars made of resistant pressed glass. The available models present different capacities for standard 76x26 mm microscopic slides, speeding up the long series of samples to be stained. Furthermore, each jar is equipped with a glass lid that avoids evaporation and dust deposition on the solution to be used.

Complete staining cuvettes

Code BPL010



Consisting of a staining tray, staining basket, and handle.

- Staining dish with pressed glass lid Code: BPL007
- Basket to hold 10 slides vertically or up to a maximum of 20 slides back-to-back, in pressed glass Code: BPL008
- Stainless steel handle for easy transfer of the basket Code: BPL009

Code	Int. dimensions (mm)	Out. dimensions (mm)
BPL010	85x70x105	75x60x90

Coplin staining jar

Code BPL005




- Made of pressed glass with lid
- 5-slot design; holds up to 5 slides in upright position or up to 10 slides back-to-back.

Code	For microscope slides (mm)
BPL005	76x26

Hellendahl staining jar

Code BPL001



- Made of pressed glass with lid
- Holds up to 8 slides vertically.
- Approx. dimensions with lid: 6.5x6.0x10 cm

Code	For microscope slides (mm)
BPL001	76x26



Schieferdecker staining jar

Code BPL003



- Made of pressed glass with lid
- Holds up to 10 porters horizontally.

Code	For microscope slides (mm)
BPL003	68x40x88

PP staining jar

Code BPL049



- Made of polypropylene, for up to 20 slides.
- Baskets available in PP, stainless steel, and pressed glass.
- Supplied with two adjustable lids, one blind and one with a window for the handle of the PP and stainless-steel baskets.

Code	Int. dimensions (mm)	Out. dimensions (mm)
BPL049	100x80x65	93x73x60

Borrel staining dishes

Code BPL029, BPL030



- With lid

Code	Ø (mm)	Height (mm)
BPL029	38	85
BPL030	40	80



Prism and tube for Borrel staining

Code **BPL032, BPL035**



Code	Description	Ø (mm)	Height (mm)	Width (mm)
BPL032	Prism for Borrel staining	-	65	23
BPL035	Borrel staining tube	38	85	-

Staining rack SS304

Code **BPL047**



- Made of bright electro-polished 304 stainless steel
- For manual staining of up to 20 slides.
- The bottom of the basket is open for quick immersion and draining.

Code	For microscope slides (mm)	Dimensions (mm)
BPL047	75x25, 75x38, 75x51	87x69x150

PP staining rack

Code **BPL048**



- Made of polypropylene
- For manual staining of up to 20 slides.
- The bottom of the basket is open for quick immersion and draining.

Code	For microscope slides (mm)	Dimensions (mm)
BPL048	75x25, 75x38, 75x51	87x69x150