



Density petroleum hydrometer, M-50 Series

Code **KLG060, KLG061, KLG062, KLG063, KLG064, KLG065, KLG066**

- Manufactured according to BS 718
- Reference temperature: 15°C
- Amplitude of the measuring range: 0.050 g/mL
- Scale division: 0.001 g/mL
- Total length: 270 mm
- Factory calibrated
- The instrument is supplied without certificate



Technical specifications

Code	KLG060	KLG061	KLG062	KLG063	KLG064	KLG065	KLG066
Range (g/mL)	0,650-0,700	0,700-0,750	0,750-0,800	0,800-0,850	0,850-0,900	0,900-0,950	0,950-1,000