



Pycnometers

KLQ004, KLQ008

Solid pycnometers are laboratory instruments which, based on the principle of liquid displacement, allow the density of solid materials of various shapes and sizes to be calculated.

Endo Glassware presents the following solid pycnometers made of borosilicate glass.

For solids, sands, gravels



Technical specifications:

Code	Capacity (mL)	Tolerance (\pm mm)	Outer \varnothing (mm)	Upper \varnothing (mm)	Funnel height (mm)	Total height (mm)
KLQ008	25	2.0	40	22	76	113
KLQ004	50	3.0	48	25	53	119