



## Stirring rod

Code LDR001, LDR002, LDR003, LDR004, LDR005, LDR011, LDR012, LDR013

Endo Glassware has practical stirring rods for stirring any type of solution contained in flasks, beakers, or other containers.


They are completely inert, which prevents contamination of samples by chemical reaction in the rod with the components of the shaken solutions, and have both ends burnt. They are available in various sizes.

### Helical rods

■ Solid glass

### Technical specifications

Code	Ø (mm)	Length (mm)	Lot
LDR001	5	200	20




### Straight rods

■ Solid glass

### Technical specifications

Code	Ø (mm)	Length (mm)	Lot
LDR002	5	200	50
LDR003	5-6	300	50
LDR004	6-7	400	50
LDR005	6-7	500	50



■ Borosilicate glass 3.3

### Technical specifications

Code	Ø (mm)	Length (mm)	Lot
LDR011	5-6	250	10
LDR012	7-8	250	10
LDR013	8-9	250	10

