



Filtres seringues

Réf. PVB001, PVB002, PVB006, PVB007, PVB011, PVB012, PVB016, PVB017, PVZ001, PVZ002, PVZ003, PVZ004, PVZ005, PVZ006, PVZ007, PVZ008, PVZ009, PVZ010, PVZ011, PVZ012, PVZ017, PVZ018, PVZ019, PVZ020, PVZ021, PVZ022, PVZ023, PVZ024

Widely used in laboratories as a practical, fast and reliable means of small volume filtration. Applications include, for example, filtration of samples for HPLC and GC, colloid separation, microanalysis and filtration of corrosive liquids. Sterile syringe filters are used for sterilisation of biological fluids, culture media and media additives.

- High-quality polypropylene body, with excellent chemical stability
- Ultrasonic welding ensures no leakage and resistance to high pressures
- Available with different membrane types and with diameters of 13 and 25 mm
- The sterile filters are gamma sterilised and are supplied individually packaged

Filtres seringues



Spécifications techniques

Référence	Membrane	Diamètre (mm)	Pore (μm)	Emballage
PVB011	Cellulose acetate	13	0,2	100 pcs / boîte
PVB012	Cellulose acetate	13	0,45	100 pcs / boîte
PVB016	Cellulose acetate	25	0,2	50 pcs / boîte
PVB017	Cellulose acetate	25	0,45	50 pcs / boîte



Filtres seringue non stériles

Spécifications techniques

Référence	Membrane	Diamètre (mm)	Pore (μm)	Emballage
PVB001	Acétate de cellulose	13	0,22	100 pcs / boîte
PVB002	Acétate de cellulose	13	0,45	100 pcs / boîte
PVB006	Acétate de cellulose	25	0,22	100 pcs / boîte
PVB007	Acétate de cellulose	25	0,45	100 pcs / boîte
PVZ001	Nylon	13	0,2	100 pcs / boîte
PVZ002	Nylon	13	0,45	100 pcs / boîte
PVZ003	Nylon	25	0,22	100 pcs / boîte
PVZ004	Nylon	25	0,45	100 pcs / boîte
PVZ005	PSE	13	0,2	100 pcs / boîte
PVZ006	PSE	13	0,45	100 pcs / boîte
PVZ007	PSE	25	0,2	100 pcs / boîte
PVZ008	PSE	25	0,45	100 pcs / boîte
PVZ009	PTFE	13	0,2	100 pcs / boîte
PVZ010	PTFE	13	0,45	100 pcs / boîte
PVZ011	PTFE	25	0,2	100 pcs / boîte
PVZ012	PTFE	25	0,45	100 pcs / boîte
PVZ017	MCE	13	0,22	100 pcs / boîte
PVZ018	MCE	13	0,45	100 pcs / boîte
PVZ019	MCE	25	0,22	100 pcs / boîte
PVZ020	MCE	25	0,45	100 pcs / boîte
PVZ021	PVDF	13	0,22	100 pcs / boîte
PVZ022	PVDF	13	0,45	100 pcs / boîte
PVZ023	PVDF	25	0,2	100 pcs / boîte
PVZ024	PVDF	25	0,45	100 pcs / boîte

PES : Polyéthersulfone, PTFE : Polytétrafluoroéthylène (Téflon), MCE : Esters de cellulose, PVDF : Fluorure de polyvinylidène