



Embudos de decantación Gilson, ENDO GLASSWARE (Familia PQD01)

Ref. PQD001, PQD002, PQD003, PQD004, PQD005, PQD006, PQD007, PQD008

- Fabricados en vidrio borosilicato
- Con llave de teflón
- Tapón de PP



Referencia	PQD001	PQD002	PQD003	PQD004	PQD005	PQD006	PQD007	PQD008
Capacidad (mL)	50	100	125	250	250	250	500	1000
Esmerilado	19/26	19/26	19/26	19/26	24/29	29/32	24/29	29/32