



Circular support

Code BJN008

- Made of plastic, fixed support without base and white colour.
- Capacity for 44 pipettes, accepts pipettes of different volumes.



Technical specifications

Code	BJN008
Capacity	32 holes Ø 12 mm, 12 holes Ø 17 mm

Holder for pipettes and thermometers

Code BJN009

- Colour white, made of PP.
- Capacity for 12 units up to 10 mm in diameter.



Technical specifications

Code	BJN009
Dimensions	205x150x220 mm

Pipette holder box

Code BJN010

- Cylindrical shape, made of polypropylene to hold pipettes of different lengths.
- The upper part of the case is screwed onto the lower part of the case for extra safety.
- The total available height can be adjusted by screwing or unscrewing the top to the desired level.



Technical specifications

Code	BJN010
Dimensions	360x70 mm



Rotating support

Code BJN011

- Holds up to 94 pipettes and has a central vertical axis for easy selection.
- Open bottom to protect pipette tips.
- Upper rotating disc has a grid pattern concentric with the lower disc to keep pipettes vertical.
- Liquid draining from the pipette tip is retained in the top of the fixed base by the outer rim.




Technical specifications

Code	BJN011
Capacity	94 holes Ø mm

Tray

Code BJN012

- Clutter-free storage of pipettes in laboratory drawers. For storage or for soaking pipettes before cleaning.
- Scalloped edges hold 16 small pipettes or 7 large ones.
- Central holder for 12 short pipettes in perpendicular direction.



Technical specifications

Code	BJN012
Dimensions	286x216x40 mm