

**TECHNICAL SPECIFICATIONS**

Trade name	Tape
Use	Particularly useful to protect any type of screw thread.
Material	Polietrafluoroetilen PTFE

**DIMENSIONAL SPECIFICATIONS**

Length	mm. 12,00
Width	mm. 12,00

**GENERAL FEATURES**

CE mark	No
Autoclavability	Yes
Sterilization	Chemical
Food contact suitable	Yes

**CHEMICAL RESISTANCE**

Categories of substances	Reference	Concen. %	T. 20°	T. 40°	T.° 60
Inorganic acid	Sulphuric acid	98	A	A	A
Organic acid	Benzoic acid	100	A	A	A
Alcohol	Ethanol	100	A	A	A
Aldehydes	Acetaldehyde	100	A	A	A
Inorganic alkali	Sodium hydroxide	50	A	A	A
Organic alkali	Aniline	100	A	A	A
Ketone	Acetone	100	A	A	A
Esther	Ethyl acetate	100	A	A	A
Hydrocarbons, halogenated	Dichloroethylene	100	A	A	A
Hydrocarbons, aromatic	Benzene	100	A	A	A
Hydrocarbons, linear	Hexane	100	A	A	A
Supersolvent	Tetrahydrofuran	100	A	A	A

Legenda (not applicable if neglected):

- **A: fair resistance;** exposure (30 days) to the chemical does not cause any damage.
- **B: sufficient resistance;** exposure causes damage of poor importance, which sometimes is only temporary.
- **C: poor resistance;** exposure to chemical is not allowed, or causes immediate permanent damage.

Results of testing with reference materials are to be considered as indication: in case of specific use it is recommended to carry out preliminary testing.

**PACKAGING**

Bag	Pcs. 5
-----	--------

**Quality Assurance**

Massimiliano Capitanio