



Motorized Stair Climbing Chair

Code **RZS001**

An electric stairlift is an innovative and efficient solution designed to improve the mobility and independence of people who have difficulty walking up and down stairs.

- Made of high-strength aluminium alloy
- Removable and washable seat and backrest cushion
- Pedal can be folded up for convenience when not in use
- Ramp with three adjustment positions, suitable for different floors
- Lift bar with three positions
- Only one person is required to operate the wheelchair
- Easy and comfortable to operate chair control
- Intelligent electromagnetic brake that provides safety for users

Load limit

160 kg



Code	RZS001
Dimensions (LxWxH)	133x48x157 cm (± 3.5 cm)
Foldaway size (LxWxH)	100x48x31 cm (± 3.5 cm)
Net weight	32 kg (± 2 kg)
Load limit	160 kg
Motor power	120 W
Motor voltage	20-29.4 V
Motor current	13.9 A
Power supply	110-240 V 50/60 Hz
Battery voltage	24 V