



Electronical balance 10 kg/0.1g, series 5161-LED

Code KBB012, KBB013

Precision balance model 5161-LED

The model 5161-LED precision balance contains the following features:

- High precision load cell.
- Resistant plastic housing.
- Big size stainless steel pan.
- High contrast easy to read LED display.
- Built-in rechargeable battery.
- Height adjustable feet.
- Full capacity taring.
- Fast response.
- Overload protection.
- Selectable measure units: kg, lb, g.
- Piece counting function.
- RS232 interface



Code	KBB012	KBB013
Maximum capacity	10 kg	20 kg
Tare range	10 kg	20 kg
Readability	0.1 g	0.1 g
Repeatability	± 0.1 g	± 0.1 g
Non-linearity	± 0.2 g	± 0.2 g
Corner error	± 0.2 g	± 0.2 g
Pan size	337x237 mm	
Dimensions (LxWxH)	370x370x125	
Power supply	AC 110V-220V	
Rechargeable battery	4V/4AH	
Operating temperature	5-35 °C	
Relative humidity	50-85 %	