



## Universal pH paper

Code **KDJ002, KDJ003, KDJ005, KDJ006, KDJ007, KDJ032, KDJ037**

- Available in different measuring ranges from pH 1 to 14 with an accuracy of one pH unit.
- Including colour scale for determining the pH of the solution.



KDJ003, KDJ006, KDJ032



KDJ002, KDJ005



KDJ007



KDJ037

Code	pH	Format
KDJ002	0-10	Box 200 strips
KDJ003	0-10	Roll 5 meters
KDJ005	1-14	Box 200 strips
KDJ006	1-14	Roll 5 meters
KDJ007	1-14	Box 80 strips
KDJ032*	1-14	3 Rolls 5 meters
KDJ037**	0-14	Box 100 strips

\*Spare for the code KDJ006 \*\*with 4 pads

## Methyl orange paper

Code **KDJ018**

- This is the simplest form of indicator paper and is used when an approximate pH determination is required.
- It does not have a colour scale, but the pH measurement is indicated by the colour change of the paper.



Code	Description	Colour change	Format
KDJ018	Methyl orange	pH < 3.1 (red), pH > 4.4 (yellow)	Box 200 strips