



1 | Baskets for general lab purposes such as transporting, autoclaving or cleaning material in ultrasonic baths.

2 | Made of stainless steel wire, they are solid and durable and present a

folding handle for an easy handling of the basket.

3 | The thin framework (5 mm) of the base and walls avoids the loss of even small material.

#### ▶ Round baskets



Code	Diameter	Height
BTS001	150 mm	150 mm
BTS002	200 mm	200 mm
BTS003	250 mm	250 mm

#### ▶ Square baskets



Code	Side	Height
BTS004	150x150 mm	150 mm
BTS005	200x200 mm	200 mm
BTS006	250x250 mm	250 mm