




Class A pipettes with bulb, glass

- Graduated according to ISO 835
- Blue graduation
- With safety bulb

Graduated pipette class AS with bulb, 0 below

Code **DGB001, DGB002, DGB003, DGB004, DGB005, DGB006, DGB007, DGB008, DGB009, DGB020**

- 0 positioned at the bottom
- Class AS with coloured band




Code	Capacity (mL)	Graduation (mL)
DGB001	1	0.1
DGB009	1	0.01
DGB002	2	0.1
DGB020	2	0.01
DGB008	5	0.05
DGB003	5	0.1
DGB004	10	0.1
DGB005	20	0.1
DGB006	25	0.1
DGB007	50	0.1

Graduated pipette class A with bulb, 0 above

Code **DGB010, DGB011, DGB012, DGB013, DGB014, DGB015, DGB016, DGB017, DGB018, DGB019**

- 0 positioned at the top
- Class A with coloured band



Code	Capacity (mL)	Graduation (mL)	Band colour
DGB010	1	0.1	yellow
DGB018	1	0.01	yellow
DGB019	2	0.01	green
DGB011	2	0.1	green
DGB017	5	0.05	red
DGB012	5	0.1	blue
DGB013	10	0.1	brown
DGB014	20	0.1	yellow
DGB015	25	0.1	white
DGB016	50	0.1	black