



Rods

Unthreaded rods for laboratory assemblies in different materials and lengths, with a diameter of 12 mm.



Chrome-plated steel

FCY002, FCY003, FCY004,
FCY005, FCY006

Code	Length (mm)
FCY002	60
FCY003	200
FCY004	350
FCY005	700
FCY006	1000

Aluminium

FCY008, FCY010, FCY011,
FCY012

Code	Length (mm)
FCY008	60
FCY010	350
FCY011	700
FCY012	1000

Stainless steel

FCY014, FCY015

Code	Length (mm)
FCY014	60
FCY015	200