




Protective gown for general use, especially suitable for use in the Healthcare, Hospital, Laboratory Equipment, Chemical, Aesthetics, Food, Cleaning, Painting, Industry and Agriculture sectors.

- Velcro fastening (5 velcro fasteners), shirt collar, elasticated cuffs.
- Individually bagged.
- 100% Latex-free product.

CE certified:

1. EU 2016/425: Relating to personal protective equipment. PPE category Category I.
2. EU 2017/745: On Medical Devices, Class I.
3. EU 10/2011: On plastic materials intended to come into contact with foodstuffs, studying their migration limits.



▶ Lab coat disposable

- Long-sleeved, single-use, Polipropylene (35 gr, PP) gown

Code	Dimensions	Color
RBB019	112x138 cm	White
RBB020	112x138 cm	Green
RBB021	112x138 cm	Blue

▶ Lab coat disposable soft

- Long-sleeved, single-use, Polipropylene (25 gr, PP) gown

Code	Dimensions	Color
RBB022	112x138 cm	White
RBB023	112x138 cm	Green
RBB024	112x138 cm	Blue