

# New release STERILIZATION

Optimize your  
sterilization processes

## 1. Sealers



ZZC002




ZZC003

- ✓ Manual cutting
- ✓ Stable base
- ✓ Preheating light indicator (ZZC003)
- ✓ Time adjustment dial 1-8 s (ZZC002)
- ✓ Airtight sealing

CODE	Sealing width	Maximum sealing length
ZZC002	3 mm	300 mm
ZZC003	15 mm	300 mm

## 2. Reels



- ✓ Protection
- ✓ Suitable for EO and autoclave
- ✓ Paper base and transparent cover
- ✓ Variety of sizes

CODE	Length	Width	Lot
ZCM014	200 m	75 mm	8
ZCM015	200 m	150 mm	4
ZCM016	200 m	200 mm	2
ZCM017	200 m	250 mm	2



# New release STERILIZATION

Optimize your  
sterilization processes

## 3. Self-adhesive bags

**IENDO**  
plasticware



- ✓ Single-use design
- ✓ Robust and durable material
- ✓ Large self-adhesive strip
- ✓ Suitable for EO and autoclave

CODE	Length	Width
ZCM009	165 mm	90 mm
ZCM010	260 mm	90 mm
ZCM011	260 mm	135 mm
ZCM012	360 mm	190 mm
ZCM013	410 mm	305 mm

## 4. Indicator tape

**Nahita**



- ✓ Heat-resistant adhesive
- ✓ Color change
- ✓ Ideal for autoclave sterilization
- ✓ Lead-free ink

CODE	Total length	Length	Lot
ZCM008	50 m	19 mm	8