



Cork Borer

Cork borer sharpener

Code BML001

- Hardened steel blade with metal handle
- Allows sharpening the end of each drill to allow it to slide easily through cork

Code	Total length (mm)
BML001	170



Cork borer

Code BML002, BML003, BML004, BML006, BML007, BML008, BML009, BML010, BML011


- Made of nickel-plated brass
- Includes solid rod for removing cork residues from the inside of the drill hole
- Drills of different gauges each with handle

Code	No. Gauge	Range (mm)
BML002	6	5/6/8/9/10/11
BML003	12	5/6.5/8/9/10/11/12.5/14/45/17/18.5/19
BML004	18	5/6.25/7.5/8.75/10/11.25/12.50/13.75/15/16.25/17.50/18.5/20/21.25/22.50/23.75/25/26.25



- Made of nickel-plated brass
- Does not include solid rod for removing cork residues from the inside of the drill hole

Code	Nº	Ø ext. (mm)
BML006	1	5
BML007	2	6.2
BML008	3	7.5
BML009	4	8.8
BML010	5	10
BML011	6	11.2



Stem for cork borers

Code BML005

- Solid rod with rubber cap to remove cork residues from the inside of the drill.

