



Portable meters

Code KDD011

pH-meter

- Equipped with pH combination electrode with integrated temperature probe.
- Easy-to-read 3.5" TFT colour display.
- Function with automatic data logging over time, ATC measurement, parameter setting.
- Data storage of up to 3,000 measurement sets.
- C-type side connection for data upload or data synchronisation.
- Equipped with pH standard solutions (pH 4.00/7.00/10.01, 50 mL, 50 mL).
- Protection class: IP5.



Code	KDD011
Model	DPH10
Electrode	D201T-H, non-refillable plastic
Parameters	pH / mV (ORP) / °C
pH	
Range	-2.00-20.00 pH
Resolution	0.1 / 0.01 pH
Accuracy	±0.002 pH ±1 digit
Temp. compensation	0.0-100.0 °C / 32.0-212.0 °F
Calibration	1-3 points autom. (USA/NIST/CN)
mV	
Range	-1999.9-1999.9 mV
Resolution	0.1 mV
Accuracy	±0.03% F.S ±1 digit
Temperature	
Range	0.0-100.0 °C / 32.0-212.0 °F
Accuracy	± 0.2 °C
Power supply	Rechargeable lithium battery