



Stainless steel cases for dissecting Instruments

Code BTQ001, BTQ002, BTQ003, BTQ004, BTQ005, BTQ006, BTQ008

- Material: stainless steel
- Suitable for autoclave sterilisation and ultrasonic cleaning
- Instruments are not included



Code	Length (mm)	Width (mm)	Height (mm)
BTQ001	170	70	30
BTQ002	180	80	40
BTQ003	200	100	40
BTQ004	220	120	50
BTQ005	250	120	50
BTQ006	300	180	50
BTQ008	400	200	100