




## Overhead stirrer, model SH-II-7C

### Code LLB007

- Digital model, with LCD screen.
- Large stirring capacity. Three stirring modes: clockwise rotation, counterclockwise rotation and clockwise/counterclockwise alternate rotation.
- Plastic housing and stainless steel head. Brushless motor. Motor overload protection. Supplied with stand and stirring rod.
- Two-year warranty.



#### Screen

Backlit LCD with blue background

#### Equipment Material

1/ Plastic housing  
2/ Stainless steel head

#### INCLUDES

Supplied with stand and stirring rod

## Specifications

Model	SH-II-7C
Code	LLB007
Max. stirring capacity (L)	40
Speed range (rpm)	100-2500
Speed accuracy (rpm)	± 1
Continuous operation	Yes
Timer	Yes
Max. viscosity (mPas)	15000
Screen	LCD
Stand plate size (mm)	200 x 300
Motor type	Brushless motor
Power supply	100-240V, 50/60Hz
Power (W)	80
Weight (kg)	7.5