



Laboratory stands

Very useful in chemistry or physics laboratories, stands are used to hold different equipments or instruments on experimental assemblies. In our catalogue, there are available stands composed of a rod and a base of different formats, as well as tripod stands.



A-shaped support stand: FCB001

Base in furnace painted cast-steel and rod in chromeplated steel. With rubber caps at the base.

Code	Base	Rod (L×Ø)
FCB001	200×30 mm	600×10 mm



Retort stand: FCB006, FCB007, FCB008, FCB009, FCB022, FCB023, FCB025

Base in cast-steel or plated-steel and rod in chromeplated steel.

Code	Colour	Material	Base	Rod (L×Ø)
FCB006	Black	Plated steel	100×175mm	500×10mm
FCB007	Black	Plated steel	125×225mm	500×10mm
FCB008	Black	Plated steel	150×250mm	500×10mm
FCB009	Black	Plated steel	150×300mm	500×10mm
FCB022	Blue	Reinforced p.steel	135×200mm	500×10mm
FCB023	Blue	Reinforced p.steel	175×300mm	500×10mm
FCB025	Blue	Reinforced p.steel	75×115mm	500×10mm



Retor stand base with rod: FCB024, FCB028, FCB029


Base in blue chromed meta and chromed rod.

Code	Colour	Material	Base	Rod (L×Ø)
FCB024	Blue	Chromed material	200×125mm	600×10mm
FCB028	Blue	Chromed material	160×100mm	500×10mm
FCB029	Blue	Chromed material	250×160mm	760×10mm



Laboratory stands

Very useful in chemistry or physics laboratories, stands are used to hold different equipments or instruments on experimental assemblies. In our catalogue, there are available stands composed of a rod and a base of different formats, as well as tripod stands.



Tripod support stand: FCB013, FCB010, FCB011, FCB012

Base in furnace painted cast-steel and rods in chromeplated steel. With rubber caps at the base.

Code	A	B	Rod	Colour
FCB013	85 mm	130 mm	500×10 mm	Blue
FCB010	Azul	170 mm	600×10 mm	Blue
FCB011	Azul	210 mm	700×10 mm	Blue
FCB012	175 mm	30 mm	800×10 mm	Blue



Tripod stand, round ring: FFB001, FFB002, FFB003, FFB004, FFB005

Ring and legs made of chromeplated steel. With rubber caps.

Code	Material	Diameter	Height
FFB001	Chromeplated steel	125 mm	200 mm
FFB002	Chromeplated steel	150 mm	200 mm
FFB003	Chromeplated steel	200 mm	200 mm
FFB004	Stainless steel	120 mm	200 mm
FFB005	Stainless steel	150 mm	200 mm



H shape support with 3 vertical rods: FCB031, FCB032

Cast iron base painted in blue, 500x250 mm

Consisting of:

- 3 vertical rods
- 2 horizontal rods
- 6 bossheads
- Rubber feet

Code	Vertical rods dim (Ø×L)	Horizontal rods dim (Ø×L)
FCB031	10x800 mm	10x500 mm
FCB032	12x800 mm	12x500 mm