



Short stem funnels, plastic

Code PNC001, PNC002, PNC003, PNC004

- Made in polypropilene (PP) autoclavable.
- Body of funnel es 60° angle.
- Provided wih external ribs.



Technical specifications

Code	Ø (cone)	L x Ø (stem)
PNC001	50 mm	50 x 7 mm
PNC002	75 mm	75 x 9 mm
PNC003	100 mm	95 x 9 mm
PNC004	150 mm	110 x 16 mm

Funnels for liquid transfer, plastic

Code PNM010, PNM011, PNM012, PNM013, PNM014, PNM015, PNM016, PNM017, PNM018

- High-Density Polyethylene (HDPE).
- Suitable for a wide range liquids and upstand on rim to reduce possible overflow.
- Hanging eye for storage.



Technical specifications

Code	Ø (cone)	Ø x L (stem)	Capacity
PNM010	80 mm	13 x 30 mm	100 mL
PNM011	110 mm	15 x 38 mm	300 mL
PNM012	130 mm	20 x 43 mm	450 mL
PNM013	140 mm	24 x 52 mm	750 mL
PNM014	180 mm	30 x 69 mm	1500 mL
PNM015	220 mm	35 x 81 mm	1750 mL
PNM016	260 mm	38 x 75 mm	4000 mL
PNM017	310 mm	43 x 80 mm	5500 mL
PNM018	400 mm	45 x 130 mm	10.000 mL