



## Self-sealing sterilization pouches

Code ZCM009, ZCM010, ZCM011, ZCM012, ZCM013

Self-sealing sterilisation pouches are an essential component of sterilisation management in medical, dental, veterinary and laboratory environments.

Designed in different sizes to ensure the safety and sterility of a variety of instruments, they are made of robust medical grade paper and are hermetically sealed with a self-adhesive strip. Their contents always remain visible through the transparent PET/CPP laminated film.



### Characteristics:

- Each bag is designed for single use only
- Robust and durable material
- Quick identification of contents
- Large self-sealing strip, easy to close
- Suitable for steam and EO gas sterilisation

| Code        | ZCM009     | ZCM010     | ZCM011     | ZCM012     | ZCM013     |
|-------------|------------|------------|------------|------------|------------|
| Length (mm) | 165        | 260        | 260        | 360        | 410        |
| Width (mm)  | 90         | 90         | 135        | 190        | 305        |
| Pack        | Box of 200 | Box of 200 | Box of 200 | Box of 200 | Box of 200 |