



Endo plasticware plastic tanks

Code **BTB014, BTB015, BTB016, BTB013, BTB017, BTB018, BTB019**

Thanks to the different sizes offered, these tanks can be used in all areas and are ideal for all types of storage. With reinforced edges and construction, they are strong enough to withstand heavy loads.

- Made of High-Density Polyethylene (HDPE)
- Stackable
- White colour
- Reference BTB019 has perforated walls with holes to facilitate drainage



Tall Form

Code	Capacity (L)	Int. dimensions (mm)	Ext. dimensions (mm)	Weight (g)
BTB014	20	380x280x195	440x315x200	1070
BTB015	50	560x330x280	625x365x285	2150
BTB016	80	650x380x320	725x440x330	3170
BTB013	100	535x450x430	600x500x445	4300



Low Form

Code	Capacity (L)	Int. dimensions (mm)	Ext. dimensions (mm)	Weight (g)
BTB017	10	390x290x90	430x320x95	630
BTB018	16	515x235x100	540x350x110	1070
BTB019*	16	515x325x100	540x350x110	1000

*With drainage holes

Endo plasticware plastic buckets

Code **BTB030, BTB028, BTB029**

- Made of Low-Density Polyethylene (LDPE)
- With pouring spout
- Graduated per litre
- White colour



Code	Capacity (L)	Lower Ø (mm)	Upper Ø (mm)	Height (mm)
BTB030	9	200	270	250
BTB028	12	220	302	280
BTB029	17	245	335	310