



## Axavet en-v5h vet infusion pump with heater

### Code ZMZ021

The AXAVET infusion pump with integrated heating system ensures accurate and safe delivery of fluids and medications, maximising therapeutic effectiveness and improving the patient's thermal comfort during treatment.



- Small and lightweight for easy portability
- 4.3" colour touch screen display
- Different settings available, such as: total perfusion volume, required flow rate, total perfusion time
- High sensitivity of settings (bubbles and occlusion pressure)
- PBT+PC material, corrosion resistant and easy to clean
- Long battery life up to 10 hours infusion time
- Suitable for transfusions
- Compatible with standard drippers
- IV pole holder included

## Technical specifications

<b>Code</b>	ZMZ021
<b>Cumulative volume</b>	0.1 to 9999 mL
<b>Display</b>	4.3" colour LCD
<b>Flow range</b>	0.1 – 2000 ml/h
<b>Modes</b>	1. Flow and volume 2. Total time and volume 3. Time and total flow 4. Body mass 5. Saved modes 6. Micro mode
<b>Accuracy</b>	± 5%
<b>Air bubble detection</b>	6 levels (20, 50, 100, 200, 400, 800 µl)
<b>Occlusion settings</b>	5 levels (50, 150, 300, 600, 900 mm Hg)
<b>Data review</b>	5000 records
<b>Protection class</b>	IP44
<b>Power supply</b>	100-240VAC, 50/60Hz
<b>Dimensions (L x W x H)</b>	13 x 9 x 14 cm
<b>Weight</b>	1.95 kg