



Kitasato flasks

Code PJH001, PJH002, PJH003, PJH004

- Endo Glassware's Kitasato flasks are robust, thick-walled flasks made of excellent physical and chemical resistance borosilicate glass.
- Widely used in chemistry and biology laboratories, they feature a side arm for connection to a vacuum pump or water aspirator and can be used with various types of funnels for rapid filtration of solutions.
- They also have approximate graduation markings and a large labeling surface in white.



Code	PJH001	PJH002	PJH003	PJH004
Capacity (mL)	250	500	1000	2000
Ø Base (mm)	85	100	132	165
Ø int. Neck (mm)	30	37	37	53
Height (mm)	140	190	235	255
Ø Side Arm (mm)	12	10	11	10